

海外製造品 修飾一覧

合成スケール	50 nmole~		
精製グレード	Salt-free	HPLC	PAGE

【ご注意/お問合せ】

*上記の精製タイプのいずれかからお選びいただけますが、鎖長、修飾タイプによって、精製方法の変更をお願いする場合がございます。

*海外製造品については、オンライン注文画面には参考価格が正しく表示されません。

価格はお見積希望の配列を記載のうえ、弊社DNA受注課 (oligoorder-jp@eurofins.com) へお問い合わせください。一覧に記載のなき場合もお気軽にお尋ねください。

修飾名	Designator 表記	部位
化学修飾(アミノ化,リン酸化,チオール化等)		
Acrydite	[5Acrd]	5'
Acrydite-DMT	[AcrDMT]	Internal
Amino Modifier C12	[AmC12]	5'
Amino Modifier C2 dT	[AmC2dT]	Internal
Amino Modifier C3	[AmC3]	5'
Amino Modifier C6	[AmC6]	5'
Amino Modifier C6 dT	[AmC6dT]	Internal
Amino Modifier C7	[AmC7]	3'
Phosphate	[PHO]	5', 3'
Thiol Mod C3 S-S	[ThiolC3]	3'
Thiol Modifier C6 S-S	[ThiolC6]	5'
Uni-Link AminoModifier	[AmUni]	5', Internal
蛍光修飾(FAM,HEX,TET等)		
6-Carboxyrhodamine (ester)	[6CRno]	5', 3', Internal
6-Joe, SE (ester)	[JOE]	5', 3', Internal
Alexa Fluor 350 (ester)	[Alexa350]	5', 3', Internal
Alexa Fluor 488 (ester)	[Alexa488]	5', 3', Internal
Alexa Fluor 532 (ester)	[Alexa532]	5', 3', Internal
Alexa Fluor 546 (ester)	[Alexa546]	5', 3', Internal
Alexa Fluor 555 (ester)	[Alexa555]	5', 3', Internal
Alexa Fluor 568 (ester)	[Alexa568]	5', 3', Internal
Alexa Fluor 594 (ester)	[Alexa594]	5', 3', Internal
Alexa Fluor 633 (ester)	[Alexa633]	5', 3', Internal
Alexa Fluor 647 (ester)	[Alexa647]	5', 3', Internal
Alexa Fluor 660 (ester)	[Alexa660]	5', 3', Internal
Alexa Fluor 680 (ester)	[Alexa680]	5', 3', Internal
Alexa Fluor 700 (ester)	[Alexa700]	5', 3', Internal
Alexa Fluor 750 (ester)	[Alexa750]	5', 3', Internal
Cy3.5	[CY35]	5'
Cy5.5	[CY55]	5'
Fluorescein-dT	[FLUdT]	Internal
Oregon Green 488 (ester)	[OG488]	5', 3', Internal
ROX (ester)	[ROX]	5', 3', Internal
TAMRA (STANDARD)	[TAM]	3'
TET	[TET]	5'
Tamra (ester) - 5-isomer	[TMR5]	5', 3', Internal
Texas Red-X (ester)	[TxRed]	5', 3', Internal
ダブル標識プローブ Mol. Beacon		
3' Dabcyl + 5' Cy3.5	[CY35][DAB]	5', 3'
3' Dabcyl + 5' Cy5.5	[CY55][DAB]	5', 3'
3' Dabcyl + 5' Oregon Green	[OG488][DAB]	5', 3'
3' Dabcyl + 5' ROX	[ROX][DAB]	5', 3'
3' Dabcyl + 5' TET	[TET][DAB]	5', 3'
3' Dabcyl + 5' Texas Red	[TxRed][DAB]	5', 3'

ユーロフィンジェノミクス(株) 修飾リスト(海外製造)一覧v202

修飾名	Designator 表記	部位
ビオチン化		
5'-Biotin	[Bio]	5'
Biotin PC (Photocleavable)	[BioPC]	5'
Biotin-ON	[BioON]	5', Internal
Biotin-TEG (STANDARD)	[BioTEG]	5', 3'
Biotin-dT	[BiodT]	Internal
修飾塩基		
2-Aminopurine	[2AmPu]	5', Internal
5-Bromo-dC	[5BrdC]	5', Internal
5-Bromo-dU	[5BrdU]	5', 3', Internal
5-Iodo-dU	[5IdU]	5', Internal
5-Methyl-dC	[5MedC]	5', 3', Internal
8-Oxo-dG	[8OxdG]	5', Internal
O6-phenyl-dI	[6phnIdI]	5', Internal
Pyrrolo-dC	[PyrrlodC]	5', Internal
RNA相同体(2-Oメチル)		
2'-O-Methyl-A	[2OMeA]	5', 3', Internal
2'-O-Methyl-C	[2OMeC]	5', 3', Internal
2'-O-Methyl-G	[2OMeG]	5', 3', Internal
2'-O-Methyl-U	[2OMeU]	5', 3', Internal
DNA相同体(イノシン,ウリジン等)		
3-Nitropyrrole	[NitPyr]	5', Internal
5-Nitroindole	[5NitdI]	5', 3', Internal
dK	[dK]	5', Internal
dP	[dP]	5', 3', Internal
クエンチャー		
BHQ1-dT	[BHQ1dT]	5', Internal
BHQ2-dT	[BHQ2dT]	5', Internal
Dabcyl (STANDARD)	[DAB]	3'
Dabcyl-dT	[DABdt]	Internal
QSY 7 (ester)	[QSY7]	5', 3', Internal
スペーサー		
Spacer 18	[SpC18]	5', 3', Internal
Spacer 9	[SpC9]	5', 3', Internal
Spacer C3	[SpC3]	5', 3', Internal
dSpacer	[Spd]	5', Internal
その他(S化,Dig等)		
2'3'-ddC	[23ddC]	3'
3'-dA (chain terminator)	[3dA]	3'
3'-dG (chain terminator)	[3dG]	3'
3'-dT (chain terminator)	[3dT]	3'
Cholesteryl-TEG	[CholTEG]	5', Internal
DNP-TEG	[DNP]	5', Internal
Digoxigenin (ester)	[Dig]	5', 3', Internal
Phosphorothioation	*	Internal
Puromycin	[Pmycn]	3'