


Fluor	5'	I	3'	Abs	Em
Alexa Fluor 350	•	•	•	346	442
Tide Fluor TF1	•	•	•	341	447
Oregon Green 488	•	•	•	490	514
Alexa Fluor 488	•	•	•	494	517
6-FAM	•		•	495	519
Fluorescein	•	•		495	520
Tide Fluor TF2	•	•	•	499	522
Rhodamine Green-X				504	531
TET	•			521	536
JOE	•	•	•	520	548
VIC				538	554
Alexa Fluor 532	•	•	•	530	555
HEX	•			535	556
Cy3	•			552	570
NED				546	570
Alexa Fluor 546	•	•	•	554	570
Alexa Fluor 555	•	•	•	555	572
TAMRA	•	•	•	544	576
Tide Fluor TF3	•	•	•	554	578
Rhodamine Red-X				574	594
PET				558	595
Cy3.5	•			581	596
ROX	•	•	•	575	602
Alexa Fluor 568	•	•	•	578	602
Texas Red-X	•	•	•	583	603
Tide Fluor TF4	•	•	•	588	610
Alexa Fluor 594	•	•	•	590	614
Alexa Fluor 633	•	•	•	621	639
Cy5	•			643	667
Tide Fluor TF5	•	•	•	656	670
Alexa Fluor 647	•	•	•	651	672
Cy5.5	•			675	694
Alexa Fluor 660	•	•	•	668	698
IRDye 700				689	700
Tide Fluor TF6	•	•	•	686	702
Alexa Fluor 680	•	•	•	684	707
DY-682				690	709
ATTO 700				700	719
Alexa Fluor 700	•	•	•	702	723
Tide Fluor TF7	•	•	•	756	775
Alexa Fluor 750	•	•	•	753	782
IRDye 800				767	786
DY782				783	800
Tide Fluor TF8	•	•	•	787	808



Dyes in shaded color are the new line of Tide Fluor dyes available from Eurofins MWG Operon.

Greyed-out dyes are not available from Eurofins MWG Operon. Bracketing dyes may be a suitable alternative for these hard-to-find fluors.

Tide Fluor and Tide Quencher are trademarks of AAT Bioquest, Inc. IRDye is a registered trademark of LI-COR Biosciences. Alexa Fluor is a registered trademark and NED, VIC, ROX, Texas Red, Oregon Green, Rhodamine Green, Rhodamine Red, and FAM are trademarks of Life Technologies. Cy is a registered trademark of GE Healthcare.

ユーロフィンジェノミクス株式会社

お問い合わせ (営業部) :

sales-jp@eurofins.com TEL03-6701-8080

ご注文 (カスタム DNA 受注課) :

oligoorder-jp@eurofins.com